

Zemits DuoCratus — Marissa Cheat Sheet

808 nm diode hair removal + IPL filters. Starting guidance only — confirm with Medical Director (Dr. Mukesh Arora, MD). Always test-spot. Device presets win if they differ.

FOR MARISSA
OPERATOR SHEET
JUL 2026

Golden rules: Start low → test spot → wait → treat. Endpoint = mild perifollicular redness / slight edema / warmth. No change → raise fluence **1–2 J/cm²** and reassess. Pain / whitening / blister → **STOP**, lower energy, max cooling, escalate. Darker skin / thicker hair = lower energy, longer pulse, max cooling. Never treat over tattoos, moles, or pigmented lesions (stay ≥ 1.5 cm away). White / gray / very light hair = poor diode response — do not force high energy.

1 · 808 NM DIODE — LASER HAIR REMOVAL (SAFE VS AGGRESSIVE)

Device range: fluence 1–120 J/cm² · pulse 5–400 ms · freq 1–10 Hz · cooling HIGH (head ~4°C before Ready). Area + Fitzpatrick on screen auto-presets — use table to sanity-check.

FITZPATRICK	MENU CUE	SAFE START (FIRST VISIT / TEST)	AGGRESSIVE* (AFTER PROVEN TOLERANCE)	PULSE / COOLING
I–II	White / fair	~16–22 J/cm ²	~24–30 J/cm ²	10–30 ms · Cooling HIGH
III–IV	Medium	~12–18 J/cm ²	~20–24 J/cm ²	20–40 ms · Cooling HIGH
V (VI = MD only)	Dark	~10–12 J/cm ²	~14–16 J/cm ²	30–100 ms · Cooling MAX

*Aggressive = high end of HG conservative band for coarse/dense dark hair after safe response — not a max-out. Coarse back/legs tolerate higher; face/neck stay lower. Overall diode working band often 10–60 J/cm² clinically; never jump >2 J/cm² between passes.

2 · DIODE TREATMENT FLOW (EVERY VISIT)

1. Consent + screen + Fitzpatrick + hair color/type + photos
2. Shave only (no wax/pluck 4 wks). Clean dry skin. Gel 2–3 mm for test/treat
3. Goggles on everyone · room warning · standby when paused
4. Select **Laser HR** → area → skin type → cooling HIGH
5. **Test spot** discreet area · wait 5–10 min · observe
6. Treat: tip flush, mild pressure, with hair growth · grid · ≤30% overlap · 1 pass unless MD says
7. Wipe tip with 70% ethanol every 30–40 pulses
8. Chart: area, Fitz, fluence, pulse, Hz, passes, response
9. Cold gauze + soothing gel · aftercare · standby → wait 5 min → power down

3 · SHOT COUNT TARGETS

AREA	SHOTS
Armpits	200–300 (100–150 each)
Bikini	450–650
Full bikini	600–800
Upper lip	60–100
Chin	50–250
Forearms	900–1,400
Full arms	1,200–2,000
Calves	~3,000
Full legs	5,000–6,000

Intervals: ~4 weeks diode HR. Series usually 6–8 (sometimes 10). Treat only when hair is growing.

4 · IPL FILTERS — PICK THE RIGHT CUT-OFF (DUOCRATUS)

IPL energy density typically 1–50 J/cm². Filters: 480 · 530 · 590 · 640 · 690. Use cooling gel. Test spot every new filter / new area.

FILTER	USE FOR	SAFE START	AGGRESSIVE*	NOTES
480 nm	Acne / active pustules	10–12 J/cm ²	13–15 J/cm ²	Porphyrins · Fitz I–III preferred
530 nm	Sun spots, freckles, facial redness / small vessels	12–15 J/cm ²	16–20 J/cm ²	Superficial pigment + hemoglobin · not melasma / not leg veins
590 nm	Hyperpigmentation + photorejuvenation · rosacea / telangiectasia	14–18 J/cm ²	19–24 J/cm ²	Hemoglobin & pigment · Fitz I–IV · pulse ~20–50 ms
640 nm	Hair removal (light–medium skin) · some vascular	18–22 J/cm ²	23–28 J/cm ²	Deeper melanin · Fitz I–IV
690 nm	Rejuvenation depth / coarse dark hair · dense body hair	20–26 J/cm ²	27–35 J/cm ²	Deep melanin · pulse often longer (~30–60 ms)

Marissa quick map: **590 = pigment / rejuvenation tone** · **690 = deeper / coarse hair & deeper rejuvenation work**. Prefer diode 808 for primary laser hair packages; use IPL 640/690 when protocol / MD directs.

5 · PRE-CARE (TELL CLIENT)

- No sun / tanning / self-tanner 4–8 weeks (longer if Fitz IV+)
- Shave 12–24 hrs before — **no wax/pluck** 4 weeks
- Stop retinoids / acids / bleaching agents on area ~14 days (screen)
- No Roaccutane / oral retinoids prior 6 months without MD clearance
- Clean skin — no lotion, makeup, deodorant

6 · AFTER-CARE

- Cool compress + soothing gel for redness/warmth
- SPF 50+ · avoid sun on area ~1 month · no tanning
- No hot showers, sauna, pool, heavy sweat 24–48 hrs
- Don't pick — shed hairs fall out over 1–2 weeks
- Makeup only when skin is intact / not fragile

Hard stops / call Medical Director or Ryan FNP before treating: pregnancy / breastfeeding · pacemaker · metal/implants in area · active cancer · active herpes in area · open wounds · recent Accutane · fresh tan · tattoo in field · Fitz VI without written MD parameters · blistering on test spot · any concern outside training. Absolute vs relative contraindication list lives in the full Zemits diode manual — when in doubt, don't treat.

